

TECHNICAL DATASHEET

Revision: 27.09.2024 page 1

Endo Phaco Tips Reusable MD0105.126

General Information

Intended use: The product is intended to be used during surgical procedures in the posterior eye segment for:

- emulsifying lens
- aspirating lens fragments and fluids

For use with device: Oertli phaco hand piece VG800011, Oertli surgery unit

Notified Body: DQS Medizinprodukte GmbH, 60433 Frankfurt am Main, Germany

Certification: CE0297

Classification: Class IIa, according to rule 6.1

GMDN-Code: 34901

Manufacturer: Oertli Instrumente AG, 9442 Berneck, Switzerland

Technical Data

Type of use: reusable

Processing instructions: see TN999042, instructions for processing

Sterilization cycles: see TN999042, instructions for processing

Instructions for use: TN999022

Packaging: See at the section "List of content" for the specific content of each article of this product group

Packaging size: 1 plastic bag contains 1 Phaco tip inside a phaco key

Material:

Part	Direct patient contact	In-direct patient contact	No patient contact
Phaco tip	Titanium		
Phaco key			Polysulfone (PSU)

Table 1: material specification

PVC: No

Phthalates: No

Animal origin: No

DESCRIPTION OF MODEL

Revision: 27.09.2024 page 2

List of content:

Article Number	VV800050
Description	Endo Phaco tip 19G, 15°, reusable
Packaging:	Inner: Phaco key Outer: plastic bag
Dimensions	Phaco tip: length: 31 mm / Ø tip: 1.1 mm / angle: 15° Phaco key: length: 30 mm / Ø: 12 mm



Table 2: list of content

Article Number	VV800051
Description	Endo Phaco tip 20G, 15°, reusable
Packaging:	Inner: Phaco key Outer: plastic bag
Dimensions	Phaco tip: length: 31 mm / Ø tip: 0.9 mm / angle: 15° Phaco key: length: 30 mm / Ø: 12 mm



Table 3: list of content