

TECHNICAL DATASHEET

Revision: 03.07.2024 page 1

Diathermy Tips Reusable MD0105.122

General Information

Intended use: The product is intended to be used during surgical procedures in ophthalmic surgery for:
- coagulation and coaptation of tissues and vessels

For use with device: Oertli HF handpieces VE201710 or VE201712 and Oertli surgical platform

Notified Body: DQS Medizinprodukte GmbH, 60433 Frankfurt am Main, Germany

Certification: CE0297

Classification: Class IIb according to Rule 6.1: 9.2

GMDN-Code: 47848

Manufacturer: Oertli Instrumente AG, 9442 Berneck, Switzerland

Technical Data

Type of use: Reusable

Sterilization method: See TN999042

Instructions for use: TN999097

Packaging: See table 2, "List of content"

Packaging size: 1 cardboard box contains 1 diathermy tip

Material:

Part	Direct patient contact	In-direct patient contact	No patient contact
Diathermy tip canula	Stainless steel, platinum-iridium, PTFE tube		
Diathermy tip housing			Stainless steel, titanium, PMMA, beryllium-copper

Table 1: material specification

PVC: No

Phthalates: No

Animal origin: No

DESCRIPTION OF MODEL

Revision: 03.07.2024 page 2

List of content

Article Number	VE201722	
Description	Diathermy tip Eraser, 50° beveled, reusable	
Packaging	Cardboard box with foamed material for protection	
Dimensions	Length:	46 mm
	Diameter:	1.1 mm
	Angle:	50°




Table 2: List of content