

TECHNICAL DATASHEET

Revision: 03.07.2024 page 1

Diathermy Forceps Reusable MD0105.135

General Information

Intended use: The product is intended to be used during surgical procedures in ophthalmic surgery for:
- coagulation and coaptation of tissues and vessels

For use with device: Oertli HF handpieces VE201710 or VE201712 and Oertli surgical platform

Notified Body: DQS Medizinprodukte GmbH, 60433 Frankfurt am Main, Germany

Certification: CE0297

Classification: Class IIb according to Rule 6.1: 9.2

GMDN-Code: 47848

Manufacturer: Oertli Instrumente AG, 9442 Berneck, Switzerland

Technical Data

Type of use: Reusable

Sterilization method: See TN999042

Instructions for use: TN999100

Packaging: See table 2, "List of content"

Packaging size: 1 cardboard box contains 1 diathermy forceps

Material:

Part	Direct patient contact	In-direct patient contact	No patient contact
Diathermy forceps tip	Stainless steel		
Diathermy forceps arms			Stainless steel, ECTFE
Diathermy forceps housing			Stainless steel, brass, PTFE, titanium

Table 1: material specification

PVC: No

Phthalates: No

Animal origin: No

DESCRIPTION OF MODEL

Revision: 03.07.2024 page 2

List of content

Article Number	VE203902	
Description	Diathermy forceps bipolar, curved, reusable	
Packaging	Cardboard box with foamed material for protection	
Dimensions	Length:	115 mm
	Diameter:	8.7 mm



Table 2: List of content